



Clinical Color Matching Guideline

Blank kleur	EW		BW		HV		LV		BWM		HVM		LVM		3D pro-HW	
Stompkleur	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm	0.4mm	0.6mm
Bleach	BL1+	BL1+	BL1+	BL1+	BL2	BL2+	BL3+	BL3 +	BL2+	BL2+	BL3 +	BL3 +	BL4 +	BL4 +	BL2+	BL1
A1	BL2	BL2	BL3	BL2+	BL3	BL3+	BL4+	BL4+	BL3+	BL3+	BL4+	BL4	BL4	BL4	BL2	BL2
A2	BL4+	BL3+	BL4+	BL3	BL4	BL4	A1+	A1+	BL4+	BL4+	A1+	A1+	A1+	A1+	BL3	BL3
A3	BL4+	BL3	BL4	BL4+	A1+	A1+	A2+	A2+	BL4	BL4	A1	A1	A1	A1	BL4	BL4
A3.5	BL4	BL4	A1+	BL4	A1	A1	A2+	A2+	A1+	A1+	A2+	A1	A2+	A2+	A1+	A1+
A4	A1+	A1+	A1	A1+	A2+	A2+	A3+	A2	A1	A1+	A2+	A2+	A2+	A2+	A1+	A1+
Brown dark	A1	A1+	A2+	A1	A2	A2+	A3	A3+	A2+	A1	A2	A2	A2	A2	A1	A1+

Born for veneer

Brilliant area Natural area Dark cover area

Aidite[®]

Make People Healthier and More Beautiful

Aidite (Qinhuangdao) Technology Co.,Ltd.

Tel : 0086-335-8587898

Fax: 0086-335-8587198

Email: info@aidite.com